English

GALIT

EXPERTS IN DENTISTRY

DENTAL UNITS



Expertise and innovation for modern dentistry

GALIT – EXPERTS IN DENTISTRY

GALIT is a European manufacturer of dental and medical equipment with more than 25 years of experience. Our solutions combine advanced engineering, high-quality European components, ergonomic design, and the flexibility to meet the evolving needs of modern dental practices.

BUILT ON EXPERTISE

Our expertise covers every stage — from product development and engineering to manufacturing and final assembly. This integrated approach allows us to create reliable, durable, and flexible dental solutions tailored to the needs of healthcare professionals worldwide.

GALIT VALUES:

QUALITY — We are committed to internationally recognised quality standards, ensuring reliability and consistency in every product we manufacture.

TRUST — We build lasting relationships through transparency, integrity, and mutual respect.

CARE — We put people at the heart of what we do, creating solutions with the needs of dental professionals and patients in mind.

INNOVATION — We continuously improve our products and embrace new technologies to create smarter, more efficient dental solutions.

PARTNERSHIP — We work closely with customers, dealers, and industry partners to create long-term value and achieve shared success.

DEDICATION — We approach every project with attention to detail, responsibility, and a commitment to delivering excellent results and reliable support.







Aesthetics and functionality.

GALIT OMNIPRATIQUE is a versatile dental unit that combines ergonomics, modern technology, and elegant design. Its wide range of configuration options allows it to be tailored to different workflows and the needs of modern dental practices.

Key Features:

- High-quality European components and durable materials — polyurethane, aluminium, steel, ceramics, and a resistant powder coating.
- Up to 5 instruments on the dentist element and up to 6 instruments on the assistant element.
- Individual choice of upholstery and patient chair, a dentist element with whip-arm or hanging-hose instrument delivery, and a membrane or touchscreen control panel.
- Swivelling ceramic spittoon for optimal access and added comfort.
- Flexible configuration options for ergonomic work and a high level of patient comfort.
- Service-friendly design with easy access to components during maintenance.







Perfection in every detail.

GALIT PRO is a premium dental unit that combines refined design with high functionality, well-considered ergonomics, and modern technology. Flexible configuration options allow it to be tailored to different workflows and modern clinic concepts.

Key Features:

- Elegant design created in collaboration with an Italian designer.
- High-quality European components for reliable performance and a long service life.
- Individual choice of upholstery, programmable chair positions, dentist element, and control panel.
- Robust construction for mounting additional equipment — a microscope, X-ray unit, or monitor.
- Two independent suction lines — high-volume evacuator and saliva ejector — for convenient and efficient work.
- Service-friendly design with easy access to components during maintenance.







Freedom from utility connections.

GALIT AUTONOME is an autonomous ALL-IN-ONE dental unit that requires no fixed water, compressed air, suction, or drainage connections. Only a power connection is needed — flexible, independent, and fast.

Key Features:

- Autonomous operation without external water, compressed air, suction, or drainage connections — only a power connection is required.
- Integrated supply system: compressor, stainless steel air tank, and powerful suction.
- Built-in H14 HEPA filter and separate containers for clean water and waste fluid.
- Flexible expansion options: microscope, X-ray unit, monitor, or air dryer.
- High-quality European components for reliable performance and a long service life.
- Individual choice of upholstery, programmable chair positions, dentist element, and control panel.
- Low-noise operation — ≤ 59 dB.

PLUG & PLAY
ONLY A POWER CONNECTION REQUIRED

LOW NOISE LEVEL
≤ 59 dB







For right- and left-handed dentists.


GALIT AMBIDEXTRE is a versatile dental unit for right- and left-handed dentists. Its well-designed rotating mechanism allows quick and easy adaptation to the dentist's preferred working position — ideal for clinics where several dentists work.

Key Features:

- Quick switching between right- and left-handed positions thanks to the flexible rotating mechanism.
- Dentist and assistant elements can be easily moved to the required side.
- Various mounting options for the operating light and dentist element — on the ceiling, floor, or patient chair, depending on the configuration.
- Dentist element with up to 5 instruments; membrane or touchscreen control panel.
- Ergonomic patient chair with individual choice of upholstery and programmable working positions.







More space. More comfort.

GALIT CONSOLE is a modern dental unit with well-considered ergonomics, high functionality, and a compact design. The direct connection between the patient chair and the water unit frees up legroom and ensures a comfortable working position during treatment.

Key Features:

- Direct connection between the patient chair and the water unit — more legroom and freedom of movement.
- Modern dentist element with up to 5 instruments.
- Membrane or touchscreen control panel.
- Ergonomic patient chair for a high level of comfort during treatment.
- High-quality European components for reliable performance and a long service life.
- Swivelling ceramic spittoon for convenient access and easy positioning.

Dentist element with whip-arm or hanging-hose delivery

Whip-arm delivery (Continental)



Hanging-hose delivery (International)



All GALIT stationary dental units are available in both versions.

Floor-standing or chair-mounted water unit

Chair-mounted water unit



Floor-standing water unit



All GALIT stationary dental units, except GALIT AUTONOME, are available in both versions. The chair-mounted water unit is attached directly to the patient chair, while the floor-standing version is installed on the floor next to the chair.



GALIT NO STRESS (CART AUTONOME and MOBI SUC EASY)





More freedom. Fewer connections.

GALIT NO STRESS is a mobile solution for flexible treatment set-ups without fixed water, compressed air, suction, or drainage connections. The system combines the CART AUTONOME mobile dental unit with a mobile suction unit — with or without a spittoon. Only a power connection is required.

Key Features:

- No fixed water, compressed air, suction, or drainage connections required.
- Combines CART AUTONOME with the MOBI SUC EASY or MOBI SUC COMBI mobile suction unit — without or with a spittoon.
- Mobile equipment that is easy to move and use in different treatment rooms.
- Built-in H14 HEPA filter for hygienic operation.
- No utility lines in the working area — more freedom of movement and a tidy workplace.
- Low noise level — ≤ 54 dB (CART AUTONOME) and ≤ 58 dB (MOBI SUC EASY).

PLUG & PLAY
ONLY A POWER CONNECTION REQUIRED

LOW NOISE LEVEL
≤ 58 dB



GALIT CART-5 and MOBI SUC COMBI





Compact. Mobile. Ergonomic.

GALIT CART-5 is a compact mobile dental unit designed for flexibility, efficiency, and ergonomic workflow. Its space-saving design and easy mobility make it an ideal solution for modern dental practices.

Key Features:

- Mobile dental unit equipped with up to 5 instruments.
- Smooth height adjustment for optimal positioning and ergonomic use.
- Integrated control panel — membrane or touchscreen — for convenient control of instruments and unit functions.
- Removable, autoclavable instrument holders for easy cleaning and hygiene.
- Water bottle system with switch to central water supply.
- Compact mobile design with smooth-rolling, lockable caster wheels.
- Smooth, easy-care surfaces for quick and easy disinfection.







Where mobility meets independence.

GALIT CART AUTONOME is a mobile, autonomous dental unit with an integrated compressor and water supply system. Requiring only a power connection, it offers maximum flexibility and freedom in everyday dental practice.

Key Features:

- Built-in oil-free compressor with a stainless steel tank.
- Integrated bottle system for autonomous water supply.
- No external water or compressed air connections required.
- Low-noise operation (≤ 54 dB) for a comfortable treatment environment.
- Membrane or touchscreen control panel.
- Removable, autoclavable instrument holders for easy cleaning and hygiene.
- Easy mobility with smooth-rolling, lockable caster wheels.
- Designed for easy cleaning and disinfection.

PLUG & PLAY
ONLY A POWER CONNECTION REQUIRED

LOW NOISE LEVEL
≤ 54 dB







ALL-IN-ONE. Mobile and convenient.

GALIT CART ALL-IN-ONE is a mobile, autonomous dental unit with an integrated compressor, water supply, and suction system. Requiring only a power connection, it provides a complete, ready-to-use solution for maximum flexibility in any dental workspace.

Key Features:

- Built-in oil-free compressor with a stainless steel tank.
- Integrated suction system with HEPA filtration.
- Integrated bottle system and 5-litre waste container for autonomous operation.
- Low-noise operation (≤ 59 dB) for a comfortable treatment environment.
- Compact, space-saving design.
- Membrane control panel.
- Up to 6 individually configurable instruments.
- Easy mobility with smooth-rolling, lockable caster wheels.
- Designed for easy cleaning and disinfection.

PLUG & PLAY
ONLY A POWER CONNECTION REQUIRED

LOW NOISE LEVEL
≤ 59 dB



GALIT CART-S and CART-5





Compact. Mobile. Powerful.

GALIT CART-S is a mobile suction unit designed for dental surgery and implantology. It provides reliable suction for efficient removal of fluids and aerosols, helping maintain a clean and dry operating field throughout the procedure. Its compact, mobile design allows flexible positioning in any dental workspace.

Key Features:

- Stable suction performance for surgical and implantology procedures.
- Clean and dry operating field for maximum precision.
- Compact, space-saving design.
- Smooth-rolling caster wheels for easy positioning during treatment.
- Medical-grade housing designed for easy cleaning and disinfection.



MOBI SUC EASY and MOBI SUC COMBI





Mobile. Powerful. Autonomous.

MOBI SUC EASY and MOBI SUC COMBI are autonomous mobile suction units designed for efficient removal of fluids and aerosols during dental procedures.

MOBI SUC EASY comes without a spittoon, while MOBI SUC COMBI features an integrated ceramic spittoon.

Key Features:

- High suction capacity of up to 300 L/min.
- Low-noise operation (≤ 58 dB) for a comfortable treatment environment.
- Mobile frame with antistatic caster wheels and integrated brakes for easy movement and secure positioning.
- Built-in H14 HEPA filter for effective removal of infectious microparticles from exhaust air.

PLUG & PLAY
ONLY A POWER CONNECTION REQUIRED

LOW NOISE LEVEL
≤ 58 dB


Foot control options



WIRELESS FOOT CONTROL
Wireless foot control with Bluetooth connectivity and a high-performance battery, ensuring a stable connection and an operating time of up to 6 months.



MULTIFUNCTION WIRED FOOT CONTROL
Precise, infinitely variable control via pressure and slide operation.



STANDARD WIRED FOOT CONTROL
Simple and reliable operation via push button.

Control panel



MEMBRANE CONTROL PANEL
Provides intuitive control of dental instruments and unit functions, with easy adjustment of key operating parameters.



TOUCHSCREEN CONTROL PANEL
An intuitive interface, advanced functionality, and real-time parameter display provide precise and efficient control of instruments and unit functions.

Dental operating lights



B75
Powerful light for surgery and implantology. Up to 160,000 lux, Composave mode (prevents premature curing of composites).

MAIA LED
Even, shadow-free light. 5000 K, up to 45,000 lux.





EVA LED
4 colour temperatures and Composave mode (prevents premature curing of composites) for aesthetic restorations. Up to 50,000 lux.

GALIT dental units can be equipped with an operating light of your choice — for everyday treatment, aesthetic restorations, or surgical procedures.

Impressive quality. Inspiring innovation.



From concept to final assembly, the GALIT team carefully manages every stage of the manufacturing process.

For customised solutions, each project begins with dedicated technical planning and engineering to ensure precision, reliability, and consistent quality.



GALIT operates 5,000 m² of modern production facilities.

Advanced manufacturing capabilities ensure high precision, consistent product quality, and efficient order fulfilment.

With an experienced engineering and design team, state-of-the-art manufacturing equipment, and partnerships with more than 30 renowned dental brands, GALIT delivers tailored solutions for even the most demanding projects.

We build our dental units to last — just like the relationships we build with our customers and partners.

GALIT – EXPERTS IN DENTISTRY





GALIT
EXPERTS IN DENTISTRY

Manufacturer

PC Galit
6E, 15 Kvitnya Str., Baikivtsi, Ternopil region,
47711, Ukraine
+380-800-502-998

office@galit.te.ua
service@galit.te.ua

EU Authorised Representative
Galdent CZ s.r.o.
K. Vapence 448, 41503 Teplice, Czech Republic
+420-775-471-378

export@galit.te.ua
galdentcz@gmail.com

CE 2797

Version: 20261022

GALIT reserves the right to make changes to the design, construction, components and technical specifications of its products as part of continuous product development and improvement, without prior notice. Images and illustrations in this catalogue are for reference only and may show previous product versions or differ from the current product configuration.